



# Disposable Hemostatic Clips

Easy-Effective-Exactly

A Xinwell Member Company

# AteTec

## Intended use

This hemostatic clips are compatible with endoscopes ,clamps are placed in the upper and lower gastrointestinal tracts to stop bleeding caused by various reasons and wounds after endoscopic surgery.

## Products Features

- Ergonomically designed handle is more accurate and convenient to operate.
- Smooth coating design to avoid damage to the instrument channel of the endoscope .
- The handle can rotate synchronously with the clamp device, which can meet the needs of different lesions on the release angle and direction of the clamp device in surgery.
- Variety spans of 9mm、11mm、13mm and 16mm,which fully meet the different requirements. Excellent clamping strength,safety and reliable.
- The clip of the clip device adopts the material with excellent toughness, which can reduce the injury of the affected part during the operation. There will be a prompt sound after release, which can improve the safety performance of the instrument and reduce the risk in the operation.
- Clips Angle 135° and re-open /close functionality to allow precise placement and control.

Model	Diameter of Outer Tube (mm)	Clip Opening (mm)	Opening angle	Working length (mm)	Characteristic
*ATE-ZXJ-A-2.3×1950×90×9	2.3	9	90°	1950	Without Infusion Port
*ATE-ZXJ-A-2.3×2350×90×9	2.3	9	90°	2350	Without Infusion Port
*ATE-ZXJ-A-2.3×1950×90×11	2.3	11	90°	1950	Without Infusion Port
*ATE-ZXJ-A-2.3×2350×90×11	2.3	11	90°	2350	Without Infusion Port
*ATE-ZXJ-A-2.3×1950×90×13	2.3	13	90°	1950	Without Infusion Port
*ATE-ZXJ-A-2.3×2350×90×13	2.3	13	90°	2350	Without Infusion Port
*ATE-ZXJ-A-2.3×1950×90×16	2.3	16	90°	1950	Without Infusion Port
ATE-ZXJ-A-2.3×2350×90×16	2.3	16	90°	2350	Without Infusion Port

\*coming soon

**JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD**

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone  
213161 Jiangsu China

E-mail: [atemedical@xinhaigroup.com](mailto:atemedical@xinhaigroup.com) Web: [www.ate-medical.com](http://www.ate-medical.com)

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

**AteTec**